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
Rührwerk R100SD

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4.1	INTENDED USE.....	26
4.2	THE FRONT PANEL CONTROLS	27
4.3	DESCRIPTION OF THE FRONT PANEL FUNCTIONS	27
4.4	DESCRIPTION OF THE REAR PANEL ON THE R100SD	28
4.5	STIRRING OPERATION	29
5	SAFETY FUNCTIONS	29
5.1	ROTOR STUCK DETECTION	29
5.2	OVERLOAD PROTECTION.....	29
5.3	OVERTEMPERATURE PROTECTION	30
6	CLEANING AND MAINTENANCE.....	30
7	DISMANTLING, TRANSPORT AND STORAGE	30
7.1	DISMANTLING.....	30
7.2	TRANSPORT AND STORAGE.....	31
8	DISPOSAL.....	31
9	WARRANTY AND LIABILITY.....	31
10	TECHNICAL DATA	32
11	REPAIRS	32
11.1	REPAIR RETURN FORM	33
12	DECLARATION OF CONFORMITY	34

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
Instruction Manual

Stirring motor R100SD

R100SD - microprocessor controlled overhead stirrer

- high power output
- back-lit alphanumeric LE-Display
- Fuzzylogic control
- Soft motor start/stop
- Overtemperature protection

1 General Information

The CAT overhead stirrer R100SD is designed in accordance with Safety Class 1, and built and tested in accordance with DIN EN 61010.

According to these regulations, the unit is designed to meet the requirements for safe and correct operations. To maintain the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

1.1 Unpacking the instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred in transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

1.2 Scope of Delivery

Please check that the package contains the following:

1 Overhead Stirrer R100SD	Part No. 60213-00
2 Support Rods	Part No. 10402-64
1 Instruction Manual	

Accessories and stirring tools (optional):

Special Clamp to fix the unit to a stand	Part No. 60495-00
Stirring shaft protector to be fixed to the stand	Part No. 20
U-Stand	Part No. 60494-00

Please see price list for further accessories such as suitable stirring paddles

Important Note:

After reading and understanding the instruction manual you may now start operating the unit.

Important Note:

Read this instruction manual carefully before operating the instrument. Should there be any additional questions, after reading these instructions, concerning the set-up, operation or warranty, please contact either your distributor, or the manufacturer at the following address

CAT

Ingenieurbüro CAT
M. Zipperer GmbH
Etzenbach 16
D-79219 Staufen
+49-7636-7803-0

Important Note:

Please store the instruction manual in a place easily accessible to every user.

A EURO AC plug (DIN 49441 CEE 7/VII10/ 16 A 250 V) is standard on all instruments. For North America instruments feature standard US plugs (NEMA Pub.No.WDI1961 ASA C 73.1. 1961 page 8 15A 125V).

When operating the instruments in countries with different AC plug systems use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation.

The instrument is earthed as supplied. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!

2 Safety Information

2.1 Explanation of the safety information in the instruction manual

This symbol emphasizes the following notices, which have to be strictly observed by the user! Disregarding of warnings may result in impairment of serviceability as well as impairment of the user.

This symbol emphasizes prohibitions, which have to be strictly observed. Any disregard of the prohibitions may result in severe impairment of serviceability as well as severe impairment of the user.

This symbol emphasizes notices, which have to be observed by the user to assure a safe operation of the unit.

CAUTION:
Hazard of fire or explosion!

Advise for repair / maintenance

Advise for connection to the mains

2.2 Explanation of the safety information on the device

CAUTION:


This symbol indicates to read the instruction manual carefully prior to operation of the instrument. Please mark points which require special attention in your field of application so they are not overlooked. Disregarding of warnings may result in impairment of serviceability as well as impairment of the user.

2.3 General Safety Information

	Attention: Please comply with all safety and accident-prevention regulations applicable to laboratory work.
	Attention: Use great caution when working in the vicinity of highly flammable or explosive substances. The instrument is not explosion-proof.
	Attention: Only instructed users may operate the instrument.
	Attention: When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.
	Turn the power switch off when not using the instrument and before disconnecting it from its AC power outlet.
	Important Note: Make sure that all safety information of the instrument respectively on the device itself are clearly visible during operation
	Attention: Do not open the instrument. Repairs are only to be carried out by trained service technicians.
	Always connect the instrument to an earthed AC power outlet.
	Ensure that the unit is firmly attached to a stand.

3 Set-up


3.1 Mounting of Support Rods

Screw the two support rods (A) into the threaded holes at the back of the overhead stirrer.

Use the chuck (B) to firmly tighten them.


3.2 Attach the Overhead Stirrer to a Stand

Mount the special clamp to the rod (C) of the stand (D).

Insert the support rods (A) of the over head stirrer into the openings of the special clamp (C) and turn the knurled screws (E) until the drive unit is firmly attached to the stand.

3.3 Attachment of Stirring Paddle

Stirring paddle (G)

Insert the stirring paddle (G) into the chuck (F) (as shown in the picture) and firmly tighten it in the desired position. Use the chuck key (B) to do so. The immersion depths of the stirring tool may be adjusted by loosening the chuck (F) and axially moving the stirring tool.

chuck (F)


Note:

To avoid damage of the chuck we recommend inserting the stirring rod minimum 50 mm into the hollow shaft of the stirring motor. For general use we recommend stirring rods with a length of 500 mm.

To achieve a balanced run of the stirring paddle please evenly fasten the chuck (F) via the 3 bores in the chuck by using the chuck key (B).

3.4 Assembling the Shaft Protector

Use a stirring shaft protector to provide protection against injury when working with the unit.

Attach the fastening rod (H) of the shaft protector (M) in front of the chuck to the stand; e.g. use cross over clamp part no. 60492-00) to fix the stirring shaft protector. Now bring the pockets (K) to the screws (J) and turn the stirring shaft protector to the left to the limit stop. The screw (L) is used to adjust the length of the stirring shaft protector.

Check that the stirring shaft protector is held in position securely prior to each use and also at regular intervals. The position of the stirring shaft protector must only be adjusted when the unit is stationary and the power supply is disconnected.

4 Set-up and Operation


4.1 Intended Use

The R100SD is designed for use in chemical and biological laboratories of industrial enterprises, universities and pharmacies. It is suitable for stirring and mixing of liquids from medium to high viscosity. For correct use the device must be fixed to an appropriate stand

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.

4.2 The Front Panel Controls

The innovative key-pad makes data entry easy and incorporates a 2-line alpha-numeric LE-display (1) which shows all relevant data. The entry and changing of the motor speed is set with the speed control knob (3).


(view on the R100SD front panel)

4.3 Description of the front panel functions

LCD display (1)	<i>The LCD display shows all relevant system data</i>
Power LED (4)	<i>indicates that unit is powered up</i>
Load LED (2)	<i>indicates when the unit is not able to hold set motor speed due to too high torque taken from motor shaft. Under this condition the unit will reduce the real motor speed so that the maximum allowed motor torque is not exceeded.</i>
Speed control knob (3)	<i>is used to set the desired motor speed</i>
Power on/off button(5)	<i>Press this key to turn the unit on or off.</i>

4.4 Description of the rear panel on the R100SD

(rear view of the R100SD)

On the rear panel of the R100SD, you will find 2 threads (a) to attach the support rods (A) and the mains power supply cable (12).

Make sure that the voltage printed on the R100SD corresponds to the voltage from your mains !

4.5 Stirring operation

CAUTION!

The units are not to be used in environments with danger of explosion.

CAUTION!

The unit is not to be used without supervision.

CAUTION!

When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

- Press the power on/off button (5)
- Set the desired motor speed using the speed setting knob (3). The digital value of the set motor speed can be read out in the top line of the display (1)

The R100SD has a soft start feature which helps to prevent the liquid from splashing.

Attention:

Avoid splashing of liquids by using a suitable vessel, lower liquid level. Make sure that suitable protective clothing and eye-wear is used.

This instruction sheet does not purport to address all of the safety problems which might result from the use of this instrument, chemicals, reagents, apparatus or equipment employed in any specific test or protocols. It is the responsibility of the user to consult their authorised safety advisors and establish appropriate health and safety practices and then determine the application of regulatory limitations prior to use.

5 Safety Functions

5.1 Rotor stuck detection

If the torque increase versus time exceeds a certain limit the unit will assume that the rotor was suddenly blocked and will switch off immediately.

Cut-off limit sensitivity can be adjusted in the Set-up menu (→ Cut off limit)

5.2 Overload protection

If the torque or power drawn from the stirring motor exceeds the allowed range, the unit will automatically reduce the motor speed until the torque is reduced to the allowed range. If this happens the LED Load is switched on.

N.B.:

As long as the Load LED is on, the motor is not able to hold the set motor speed.

5.3 Overtemperature protection

If the temperature inside the unit exceeds the maximum operating temperature, the unit will automatically reduce power output (speed and torque) in order to get back to a safe operating point (= fold back characteristic). This prevents the unit to be destroyed by constant overload operation. This overload condition is also signalled by the red 'Load' LED.

Furthermore a cooling fan on the rear side of the instrument is automatically switched on if required.

If the temperature inside the R100SD exceeds 70°C by any means, the unit will shut down immediately. This safety feature is fully "redundant".

Important Note:

Before re-starting the unit after cooling down switch it off and then on again by using the on/off button.

6 Cleaning and Maintenance


The outer casing consists of painted or anodised aluminium and a chemically resistant splash-proof membrane key-pad and therefore easily cleaned with warm water and any suitable liquid laboratory detergent. Do not use steel wool or any similar plastic wool sponge to clean the unit.

In case of malfunction do not attempt to repair the unit. There are no user-serviceable parts in this instrument.

The CAT R100SD should only be opened and repaired by authorised service personnel. Any work on the electronics inside the unit should only be carried out by knowledgeable, trained personnel. Any attempt by the user to repair the unit will immediately render the guarantee null and void. Please contact your local distributor in the event of a problem.

7 Dismantling, Transport and Storage

7.1 Dismantling

1. Switch the unit off.
2. Disconnect the unit from the mains.
3. Remove any glass beakers and any other equipment around the overhead stirrer.
4. Remove the stirring paddle (G) by loosening the chuck (F) with the chuck key (B).
5. Loosen the special clamp (C) and remove the clamp from the support rods.
6. Unscrew the support rods (A).
7. Now you may remove the instrument from the working area.

7.2 Transport and Storage

Prior to transport:

Switch the instrument off and unplug the power supply.

Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.

Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).

Do not subject the instrument to mechanical shocks or vibration during transporting it.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)

8 Disposal

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, sort according to materials: metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

9 Warranty and Liability

The manufacturer agrees to either repair, or replace, at the manufacturer's discretion, any defects in materials or workmanship which develop within 24 months of the delivery of this product to the original user. In the event of replacement, the replacement unit will be guaranteed for the remainder of the original twenty-four (24) month period or ninety (90) days, whichever is longer.

If this product should require service, contact your local distributor or manufacturer for necessary instructions.

This guarantee will not apply if the defect or malfunction was caused by accident, neglect, unreasonable use or fitness for a particular purpose, which extend beyond the description and period set forth herein.

The manufacturer's sole obligation under this guarantee is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damages of any kind, resulting from use or possession of the product.

Attention:

The user has to determine, if the instrument is suitable for his specific application. If there are any further queries, contact your local dealer or the manufacturer direct.

10 Technical Data

Type: R100SD	Specifications
Electrical power requirements	110/230V, 300W, 50/60 Hz (see R100SD rear panel)
LED-Display	bright, back-lit LED display, showing set and real values of motor speed.
Motor specification	<ul style="list-style-type: none"> ▪ microprocessor controlled for easy operation and accuracy ▪ soft-start and soft-stop logic prevents splashing of liquid ▪ precise adjustable motor speed from 40 to 2000 rpm ▪ max. 180 watts shaft power output ▪ extra long life, high performance, brush motor ▪ silent belt drive power transmission ▪ max. torque: 100 Ncm (20 min) ▪ continuous torque output (at 1100 rpm): 100 Ncm
Multistep safety system	<p>detects and protects against the following dangerous situations:</p> <ul style="list-style-type: none"> ▪ rotor stuck/block detection ▪ motor overtemperature monitoring ▪ self test of all safety functions after switching on
Permissible ambient temperature	5 - 40 °C within operating area, e.g. fume cupboard
Permissible humidity	80 % RH
Safety class acc. to DIN 40050	IP42
Case dimensions	77.5 x 190 (370 including rod) x 170 mm Length including chuck 230 mm
weight	5 Kg

11 Repairs

When returning instruments for repair that have come into contact with hazardous substances, please:

Fill in attached "Repair Return Form"

Provide precise information on the relevant medium

Take protective measures to ensure the safety of our receiving and maintenance personnel

Mark the package as appropriate for hazardous materials.

11.1 Repair Return Form

CONTACT/USER INFORMATION

Contact:	Phone No.
Fax No.	Email:
Billing:	Shipping:
Company	Company
Address	Address

INSTRUMENT INFORMATION

Model	Serial No.
Please describe all problems/malfunctions	

Operating Conditions (please fill in if applicable)		
Ambient Temp.	Humidity	Speed
Load	Volume	Viscosity
Temperature in °C	Sample Temperature	Operating Time
Sample Description*		

*NOTE: If the instrument was exposed to hazardous material, it must be decontaminated BEFORE returning it to CAT and an MSDS for hazardous material must be included with the instrument.

RETURN SHIPPING

UPS <input type="checkbox"/>	Air Parcel Post <input type="checkbox"/>	Collect** <input type="checkbox"/>	Other** <input type="checkbox"/>
------------------------------	--	------------------------------------	----------------------------------


**Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means.

PACKAGING INSTRUCTIONS TO RETURN A CAT INSTRUMENT FOR REPAIR

- ✓ Remove all accessories (e.g. homogniser tools, stirring paddles) from the instrument
- ✓ Clean excess testing material off the instrument/accessory
- ✓ Include MSDS sheets for all hazardous materials used with this instrument
- ✓ Pack the instrument in its original box. If the box is not available, take care to wrap the instrument and accessories with enough material to support them.
- ✓ DO NOT send pedestal stand unless there is a problem with the upright rod, clamp or base. If there is a problem with the stand remove the upright rod from the base and individually wrap each item to avoid contact with the instrument. (Applicable for overhead stirrers and homogenisers.)
- ✓ Pack the instrument and related items in a strong box for shipping. Mark the outside of the box with handling instructions.

Example: "Handle with care" or "Fragile- Delicate Instrument" and send to:

Ingenieurbüro CAT, M. Zipperer GmbH, Service Department, Etzenbach 16, D-79219 Staufen

Ingenieurbüro
M. Zipperer GmbH

Etzenbach 16,D-79219 Staufen
Tel.: +49 (0) 7636 78030
Fax: +49 (0) 7636 780345
E-Mail: info@cat-ing.de
Internet: <http://www.cat-ing.de>

12 Declaration of conformity

We declare under our sole responsibility that this product (see table) corresponds to the regulations and conforms with the standards or standardized documents

2006/95/EC

DIN EN 61010-1; VDE 0411-1:2002-08;Amendment 1; VDE 0411-1 Amendment 1:2002-11; Amendment 2;
VDE 0411-1 Amendment 2:2004-01
DIN EN 61010-2-051; VDE 0411-2-051:2004-07

2004/108/EC

DIN EN 61326-1; VDE 0843-20-1:2006-10
DIN EN 61000-6-2; VDE 0839-6-2:2006-03
DIN EN 61000-6-4; VDE 0839-6-4:2007-09
DIN EN 61000-4-2; VDE 0847-4-2:2009-12
DIN EN 61000-4-3; VDE 0847-4-3:2008-06
DIN EN 61000-4-4; VDE 0847-4-4:2005-07
DIN EN 61000-4-5; VDE 0847-4-5:2007-06
DIN EN 61000-4-6; VDE 0847-4-6:2009-12
DIN EN 55011; VDE 0875-11:2010-05
DIN EN 55014-1; VDE 0875-14-1:2010-02
DIN EN 55014-2; VDE 0875-14-2:2009-06

2006/42/EC, 93/68/EEC

DIN EN ISO 12100-1:2004-04
DIN EN ISO 12100-2:2004-04

In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.

Ingenieurbüro CAT
M. Zipperer GmbH

Manfred Zipperer, CEO

79219 Staufen, den 27. Mai 2011

Type
R100 SD 230V + 10% 50/60Hz, 300W