Operating Instructions

Stirring Motors R18, R50, R50D and R80D-PC

Inhalt

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1 General Information

The CAT overhead stirrers R18, R50, R50D and R80D-PC are designed in accordance with Safety Class 1, and built and tested in accordance with DIN EN 61010.

According to these regulations, the unit is designed to meet the requirements for safe and correct operations. To maintain the proper safety and operational functions of the instrument, the user should follow the instructions and safety guidelines in this manual.

2 Unpacking the instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred in transport be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

2.1 Scope of Delivery

Inspect the content of the package for damage or missing parts:

	1 Overhead stirrer R18,	230 V	(Part no.: 60214-00)
Or	1 Overhead stirrer R18,	115 V	(Part no.: 60214-01)
Or	1 Overhead stirrer R50,	230 V	(Part no.: 60215-00)
Or	1 Overhead stirrer R50,	115 V	(Part no.: 60215-01)
Or	1 Overhead stirrer R50D,	230 V	(Part no.: 60216-00)
Or	1 Overhead stirrer R50D,	115 V	(Part no.: 60216-01)
Or	1 Overhead stirrer R80D,	230 V	(Part no.: 60310-00)
Or	1 Overhead stirrer R80D,	115 V	(Part no.: 60310-01)
Or	1 Overhead stirrer R80D-PC,	230 V	(Part no.: 60311-00)
Or	1 Overhead stirrer R80D-PC,	115 V	(Part no.: 60311-01)
	1 Support Rod		(Part no.: 50215-02)
	1 Chuck key		(Part no.: 20)
	1 Instruction Manual		

Accessories and Stirring Paddles (optional)

Protector	(Part no.: 20618-60)
Pedestal Stand incl. rod	(Part no.: 60491-00)
Cross over clamp	(Part no.: 60492-00)
Cross over clamp made of steel	(Part no.: 60492-50)
Blade stirring paddle	(Part no.: 8B009405)
Centrifugal stirring paddle	(Part no.: 8B009211)
Propeller stirring paddle	(Part no.: 8B009105)

Optional Accessories for R80D-PC

9 Pin SUB-D cable, 1.8 m long	(Part no.: 30275-51)
USB-RS232 Adaptor	(Part no.: 30244-01)

**Important Note:**

Read this instruction manual carefully before operating the instrument. Should there be any additional questions, after reading these instructions, concerning the set-up, operation or warranty, please contact either your distributor, or the manufacturer at the following address

**Important Note:**

After reading and understanding the instruction manual you may now start operating the unit.

CAT

Ingenieurbüro CAT
M. Zipperer GmbH
Etzenbach 16
D-79219 Staufen
Tel: +49-7636-7803-0
Fax: +49-7636-7803-45
Email: info@cat-ing.de

**Important Note:**

Please store the instruction manual in a place easily accessible to every user.

A EURO AC plug (DIN 49441 CEE 7/VII10/ 16 A 250 V) is standard on all instruments. For North America instruments feature standard US plugs (NEMA Pub.No.WDI1961 ASA C 73.1. 1961 page 8 15A 125V).

When operating the instruments in countries with different AC plug systems use an approved adapter or have a qualified electrician replace the AC plug with an approved model suitable for the country of operation.

The instrument is earthed as supplied. When replacing the original AC plug, ensure that the earth conductor is connected to the new plug!

3 Safety Information

3.1 General Safety Information

Attention:

Please comply with all safety and accident-prevention regulations applicable to laboratory work.

Attention:

Use great caution when working in the vicinity of highly flammable or explosive substances. The instrument is not explosion-proof.

Attention:

Only instructed users may operate the instrument.

Attention:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

Turn the power switch off when not using the instrument and before disconnecting it from its AC power outlet.

Important Note:

Make sure that all safety information of the instrument respectively on the device itself is clearly visible during operation.

Attention:

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

Always connect the instrument to an earthed AC power outlet.

Ensure that the unit is firmly attached to a stand.

Protect the vessel in which you stir against rotating, shaking or moving.

Make sure that the unit has been turned off before exchanging the stirring elements.

Only turn the motor on and off while the stirring rod is well inserted into the media that is to be stirred.

Avoid splashing of liquids by using a tall vessel, lower level of liquid, or by using a cover on top of the container.

For protection of the user against splashing, dress accordingly and use eye protection.

Take care that the motor does not come into contact with the media to be stirred.

3.2 Explanation of the safety information in the instruction manual

This symbol emphasizes the following notices, which have to be strictly observed by the user! Disregarding of warnings may result in impairment of serviceability as well as impairment of the user.

This symbol emphasizes prohibitions, which have to be strictly observed. Any disregard of the prohibitions may result in severe impairment of serviceability as well as severe impairment of the user.

This symbol emphasizes notices, which have to be observed by the user to assure a safe operation of the unit.

CAUTION
Hazard of fire or explosion!

Advise for repair / maintenance

Advise for connection to the mains

3.3 Explanation of the safety information on the device

CAUTION:
This symbol indicates to read the instruction manual carefully prior to operation of the instrument. Please mark points which require special attention in your field of application so they are not overlooked. Disregarding of warnings may result in impairment of serviceability as well as impairment of the user.

4 Intended Use

The CAT overhead stirrers R18, R50, R50D and R80D-PC are designed for use in chemical and biological laboratories of industrial enterprises, universities and pharmacies. They are suitable for stirring and mixing of liquids from aqueous to medium viscosity. For correct use the devices must be fixed to an appropriate stand.

To ensure maximum service life, observe the specified ambient conditions (temperature and humidity) and ensure that the instrument is not exposed to a corrosive atmosphere.

4.1 General Description


The CAT laboratory overhead stirrers R18, R50 (without digital speed display), the R50D (with digital speed display), R80D (with digital display) and R80D-PC are microprocessor controlled. They are particularly suitable for all applications where solutions of average viscosity have to be stirred. The use of the instrument is advantageous where the speed adjustment of stirring processes has to be reproducible - such as when test series have to be carried out one after the other under comparable conditions, and when parallel experiments with several stirrers must be made simultaneously.

The R50, R50D and R80Dpc stirrers can also be used in the treatment of special solutions where higher viscosity values occur for short periods as the individual constituents are added.

The speed control is steplessly adjustable between 20 and 700 rpm for the R80D-PC, between 50 and 1600 rpm for stirring motors R50 and R50D; for the R18 the speed is adjustable between 110 and 2000 rpm. The electronic feedback control guarantees exact control and high precision so that fluctuations from the viscosity of the media being stirred will not occur.


5 Set-up

5.1 Mounting of Support Rod

- Screw the support rod (A) to the rear panel of the overhead stirrer with the two M5 screws (B). Use a hexagon socket screw key with AF (across flats).
- Tighten the screws thoroughly.


5.2 Attach the Overhead Stirrer to a Stand

- Mount the cross over clamp (C) to the rod of the stand.
- Insert the support rod (A) of the over head stirrer into the opening of the cross over clamp (C) and turn the knurled screw (D) until the drive unit is firmly attached to the stand.


5.3 Attachment of Stirring Paddle

R50, R50D and R80D-PC have a hollow shaft with a 10.5 mm bore. The chrome plated chuck avoids corrosion. The chuck takes any stirring rod up to 10 mm (stirring rod not included). The support rod at the back of the unit is used to attach the motor to a stand (stand not included). R18 does not have a hollow shaft but also takes any stirring rod up to 10 mm (stirring rod not included).

- Insert the stirring paddle into the chuck (E) and firmly tighten it in the desired position. Use the chuck key (G) to do so.
- To center the stirring paddle we recommend to evenly tighten the chuck via the 3 chuck bores (F).
- The immersion depths of the stirring tool may be adjusted by loosening the chuck (E) and axially moving the stirring tool.

5.4 Assembling the Shaft Protector

Use a stirring shaft protector to provide protection against injury when working with the unit.

Attach the fastening rod (H) of the shaft protector (M) in front of the chuck to the stand; e.g. use cross over clamp part no. 60492-00) to fix the stirring shaft protector. Now bring the pockets (K) to the screws (J) and turn the stirring shaft protector to the left to the limit stop.

The screw (L) is used to adjust the length of the stirring shaft protector.

Check that the stirring shaft protector is held in position securely prior to each use and also at regular intervals. The position of the stirring shaft protector must only be adjusted when the unit is stationary and the power supply is disconnected.

6 Operation

Attention:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

Caution

The unit may not be operated in explosion-prone areas.


Pull main plug as a protection against unauthorised or accidental operation whenever the work is interrupted for extended periods.


Always start with the lowest speed. Increase the speed until the speed slowly until the desired speed setting has been reached.

6.1 Description of the front panels R18 and R50


6.2 Description of the front panels R50D and R80D-PC

The R50, R50D and R80Dpc stirrers have two independent safeguards against overload:

1. Should the allowable torque be exceeded, a red warning LED will show and the electronic system automatically reduces the power. As soon as the overload is either removed or overcome the red warning LED will go off. At continuous overload and the resulting reduced power to a speed rpm < 20 the motor will be switched off.
2. Should an overload occur from overheating during continual operations or too hot of an ambient, the red LED will show and the electronic will reduce the power intake. The motor will switch off if continuously overheated. After a cool down period and switching the unit off and on the motor will operate as normal.

6.2.1 Description of the rear R80D-PC

6.3 Operation R18, R50 and R50D

Attention:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

- Insert stirring rod into chuck and secure with chuck key. To center the stirring paddle we recommend to evenly tighten the chuck via the 3 chuck bores (F).
- Insert plug into the mains outlet.
- Place stirring rod into the media to be stirred.
- Press switch to turn the device on.
- Adjust the speed by turning the speed regulator knob to the left. The speed for the R18 and R50 is reflected by the scale around the regulator knob. The speed for the R50D is indicated on the digital display. Always start with the lowest speed.
- To turn off, reduce the speed to the lowest rate, turn the switch off, and then remove the stirring rod from the media.

6.4 Operation R80D-PC

Attention:

When connecting the instrument to an AC power outlet, ensure that your local supply voltage matches that indicated on the instrument's rating plate.

- Insert stirring rod into chuck and secure with chuck key. To center the stirring paddle we recommend to evenly tightening the chuck via the 3 chuck bores (F).
- Insert plug into the mains outlet.
- Place stirring rod into the media to be stirred.
- Press switch to turn the device on.
- Adjust the speed by turning the speed regulator knob to the left. The speed for the R80D-PC is indicated on the digital display. Always start with the lowest speed.
- Connect serial cable to the device. Set the regulator knob to the left end position minimum RPM). Switch the device on to control the device over the serial port.
- To turn off, reduce the speed to the lowest rate, turn the switch off, and then remove the stirring rod from the media.


6.5 Serial Port Operation (for R80Dpc only) applicable from Version 3.1 (11.05.2007)

6.5.1 Pin out of the I/O connector

The following signals are led out on the 9-pin SUB-D connector on the rear side of the instrument:

- RS232 control signals

Pin Number	Description
1	n.c.
2	TxD (RS232 transmit data)
3	RxD (RS232 receive data)
4	n.c.
5	GND (signal ground)
6	n.c.
7	n.c.
8	n.c.
9	n.c.

Settings:

Baudrate: 1200
Data bits: 8 bit
Parity: none
Stop bits: 1

Commands:

Query of version 0,RTY,1
Change operation mode 0,WSM,1 (1 – PC, 0 – manual)
Set RPM: 0,WSE,xxx (“xxx” = 20 to 700)
Read RPM: 0,RSE,1
Write device address: 0,WSA,x (“x” = 1 to 8)

6.5.2 Programming example:**Example 1: set RPM**

1,WSM,1 ;PC mode
1,WSE,100 ;programs speed, (device address “1”)

The device replies with the following code string:

1,HS,OK ;command has been accepted and executed

Example 2: Read RPM

1,WSM,1 ;PC mode
1,RSE,1 ;Read RPM of device 1

The device replies with the following code string:

1,HS,OK,100 ;command has been accepted and executed; current RPM is 100

7 Maintenance

Attention:

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

Apart from the motor carbon brushes the instruments are maintenance-free, except from normal wear and the statistical failure rate of the components. Should the instrument require repairs, it is recommended that they be done at the CAT factory. Please send any defective instrument either to the CAT factory or your CAT agent.

Caution

Unplug the AC power cord before opening the housing.

The carbon brushes in the motor are subject to natural wear. To change the carbon brushes please carry out the following steps.

1. Remove the black cover at the side of the motor.
2. Bent up the metal clip and take out the spring.
3. Use a small screwdriver to remove the clip from the carbon and then take out the carbon.
4. Replace the carbon. Attach the spring. Bent aback the clip. Attach the black cover.

Should the tooth belt have to be changed, then:

Important Note:

The following information does not apply for overhead stirring motors R80D and R80D-PC!!!

1. Remove the hard plastic undercover.
2. Pull the tooth belt from the large cogwheel.
3. Put the new tooth belt on the small cogwheel and then onto the large cogwheel.
4. Now turn lightly the large cogwheel until the belt is in the proper position.

All bearings of the stirring motor are sufficiently lubricated so that with careful handling no maintenance is necessary. Should problems occur, please send the instrument to CAT or your CAT agent.

8 Dismantling, Transport and Storage

8.1 Dismantling

1. Switch the unit off.
2. Disconnect the unit from the mains.

- Remove any glass beakers and any other equipment around the overhead stirrer.
- Remove the stirring tool by loosening the chuck
- Loosen the cross over clamp and remove the clamp from the support rod
- Unscrew the support rod.
- Now you may remove the instrument from the working area.

8.2 Transport and Storage

Prior to transport:

Switch the instrument off and unplug from power supply.

Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.

Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).

Do not subject the instrument to mechanical shocks or vibration during transporting it.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)

9 Warranty and Liability

The manufacturer agrees to correct for the original user of this product, either by repair, or at the manufacturer's election, by replacement, any defects in material or workmanship which develop within 24 months after delivery of this product to the original user. In the event of replacement, the replacement unit will be warranted for the remainder of the original twenty-four (24) months period of ninety (90) days, whichever is longer.

If this product should require service, contact your distributor or manufacturer for necessary instructions.

This warranty shall not apply if the defect or malfunction was caused by accident, neglect, unreasonable use, improper service, or other causes not arising out of defects in material or workmanship.

There are no warranties, expressed or implied, including, but not limited to, those of merchantability or fitness for a particular purpose, which extended beyond the description and period set forth herein.

The manufacturer's sole obligation under this warranty is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damages of any kind resulting from use or possession of the product.

10 Technical Data

R80D-PC and R80D	
Motor	Permanently generated direct –current motor
Speed Range	20 – 700 1/min
Nominal Rating	120 Watt input Short Term Overload 30%
Torque over entire speed range	80 Ncm
Capacity of Chuck	0.5 – 10.5 mm
Viscosity of Medium	Average values (water to wallpaper paste)
Mains Voltage	200 -250 V, 50/60 Hz or 110 V +/- 10%, 50/60 Hz
Permissible period of operation	Continuous operation may require service intervals
Permissible ambient temperature	+ 5°C to 40°C
Diameter of support rod	13 mm
Length of support rod	196 mm
Weight	2.6 kg
Dimensions	65 x 192 (incl. 55 mm chuck) x 120 mm (w x h x d)

R50 and R50D	
Motor	Permanently generated direct –current motor
Speed Range	50 – 1600 1/min
Nominal Rating	120 Watt input Short Term Overload 30%
Torque over entire speed range	55 Ncm
Capacity of Chuck	0.5 – 10.5 mm
Viscosity of Medium	Average values (water to wallpaper paste)
Mains Voltage	200 -250 V, 50/60 Hz or 110 V +/- 10%, 50/60 Hz
Permissible period of operation	Continuous operation may require service intervals
Permissible ambient temperature	+ 5°C to 40°C
Diameter of support rod	13 mm
Length of support rod	196 mm
Weight	2.4 kg
Dimensions	65 x 192 (incl. 55 mm chuck) x 120 mm (w x h x d)

R18	
Motor	Permanently generated direct –current motor
Speed Range	110 – 2000 1/min
Nominal Rating	120 Watt input Short Term Overload 30%
Torque over entire speed range	18 Ncm
Capacity of Chuck	0.5 – 10.5 mm
Viscosity of Medium	Aqueous fluids
Mains Voltage	200 -250 V, 50/60 Hz or 110 V +/- 10%, 50/60 Hz
Permissible period of operation	Continuous operation may require service intervals
Permissible ambient temperature	+ 5°C to 40°C
Diameter of support rod	13 mm
Length of support rod	196 mm
Weight	2.2 kg
Dimensions	65 x 192 (incl. 55 mm chuck) x 120 mm (w x h x d)

11 Disposal

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, sort according to materials: metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

12 Repairs

When returning instruments for repair that have come into contact with hazardous substances, please:

Fill in attached "Repair Return Form"

Provide precise information on the relevant medium

Take protective measures to ensure the safety of our receiving and maintenance personnel

Mark the package as appropriate for hazardous materials.

12.1 Repair Return Form

Important: <ol style="list-style-type: none"> 1. Calibration of the device is done automatically. 2. You will receive an estimate for the expected repair costs. 3. Repairs will be done after release.

CONTACT/USER INFORMATION	
Contact:	Phone No.
Fax No.	Email:
Billing:	Shipping:
Company	Company
Address	Address

DEVICE INFORMATION	
Model	Serial No.
Please describe all problems/malfunctions	

Operating Conditions (please fill in if applicable)		
Ambient Temperature	Humidity	Speed
Load	Volume	Viscosity
Temperature in °C	Sample Temperature	Operating Time
Sample Description*		

*NOTE: If the device was exposed to hazardous material, it must be decontaminated BEFORE returning it to CAT and an MSDS for hazardous material must be included with the device.

RETURN SHIPPING

UPS <input type="checkbox"/>	Air Parcel Post <input type="checkbox"/>	Collect** <input type="checkbox"/>	Other** <input type="checkbox"/>
------------------------------	------------------------------------------	------------------------------------	----------------------------------

**Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means

PACKAGING INSTRUCTIONS TO RETURN A CAT DEVICE FOR REPAIR

- ✓ Remove all accessories (e.g. homogenizer tools, stirring paddles) from the device
- ✓ Clean excess testing material off the device/accessory
- ✓ Include MSDS for all hazardous materials used with this device.
- ✓ Pack the device in its original box. If the box is not available, take care to wrap the device and accessories with enough material to support them.
- ✓ Do NOT send pedestal stand unless there is a problem with the upright rod, clamp or base. If there is a problem with the stand remove the upright rod from the base and individually wrap each item to avoid contact with the device. (Applicable for overhead stirrers and homogenisers)
- ✓ Pack the device and related items in a strong box for shipping. Mark the outside of the box with handling instructions

Example: "Handle with care" or "Fragile-Delicate Instrument" and send to:

Ingenieurbüro CAT, M. Zipperer GmbH, Service Department, Etzenbach 16, D-79219 Staufen

Etzenbach 16,D-79219 Staufen
Tel.: +49 (0) 7636 78030
Fax: +49 (0) 7636 780345
E-Mail: info@cat-ing.de
Internet: <http://www.cat-ing.de>

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We declare under our sole responsibility that this product (see table) corresponds to the regulations

2006/95/EC

DIN EN 61010-1; VDE 0411-1:2002-08; Amendment 1; VDE 0411-1 Amendment 1:2002-11; Amendment 2;
VDE 0411-1 Amendment 2:2004-01
DIN EN 61010-2-051; VDE 0411-2-051:2004-07

2004/108/EC

DIN EN 61326-1; VDE 0843-20-1:2006-10
DIN EN 61000-6-2; VDE 0839-6-2:2006-03
DIN EN 61000-6-4; VDE 0839-6-4:2007-09
DIN EN 61000-4-2; VDE 0847-4-2:2009-12
DIN EN 61000-4-3; VDE 0847-4-3:2008-06
DIN EN 61000-4-4; VDE 0847-4-4:2005-07
DIN EN 61000-4-5; VDE 0847-4-5:2007-06
DIN EN 61000-4-6; VDE 0847-4-6:2009-12
DIN EN 55011; VDE 0875-11:2010-5
DIN EN 55014-1; VDE 0875-14-1:2010-02
DIN EN 55014-2; VDE 0875-14-2:2009-06

2006/42/EC

DIN EN ISO 12100-1:2004-04
DIN EN ISO 12100-2:2004-04

In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.

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M.Zipperer GmbH