**CAT**

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## *Instruction Manual*

---

*Compressed-Air Driven Homogenising System  
XP28, XP44*

## Content

1	General Information .....	5
2	Unpacking the instrument .....	5
1	.....	5
2.1	Scope of delivery .....	5
2.2	Homogenising tools and accessories (to be ordered separately): .....	5
2.3	Recommended Accerssories not available at CAT .....	6
3	Safety Information .....	8
3.1	Explanation of the safety information in the instruction manual .....	8
3.2	Explanation of the safety information on the device .....	8
3.3	General Safety Information .....	8
4	Set-up and Operation.....	9
4.1	Intended Use.....	9
4.2	Description XP44/XP28.....	9
4.3	Set-up .....	9
4.4	Air supply requirements .....	10
4.5	Operating Conditions.....	10
4.6	Assembly of support rod .....	10
4.7	Attaching the drive unit to a pedestal stand .....	10
4.8	Inserting the homogenizing tool.....	11
4.9	Operation XP44 .....	11
4.10	Operation XP28.....	12
4.11	Exchanging the homogenizing shafts .....	12
4.12	Exchanging the generator:.....	12
5	Maintenance and Cleaning.....	13
5.1	Cleaning the drive unit XP44/XP28.....	13
5.2	Cleaning the homogenising tool .....	13
5.3	Maintenance of drive unit XP44/XP28 .....	13
5.4	Maintenance homogenizing tools .....	13
6	Dismantling, Transport and Storage .....	14
6.1	Dismantling.....	14
6.2	Transport and Storage.....	14
7	Disposal.....	14
8	Warranty and Liability .....	15
9	Technical Data .....	15

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## 1 General Information

The CAT compressed- air driven homogenising system XP28/XP44 is designed in accordance with Safety Class 1 and built and tested in accordance with DIN EN 61010.

According to these regulations the unit is designed to meet the requirements for safe and correct operations. To maintain the proper safety and operational functions of the instrument the user should follow the instructions and safety guidelines in this manual.

## 2 Unpacking the instrument

Unpack the instrument carefully and check to see that it is not damaged. It is important that any damage incurred in transport to be recognized at the time of unpacking. Notify your carrier or forwarding agent immediately in case of such damage.

### 2.1 Scope of delivery

Please check that the package contains the following:

Type	Description	Part No.
XP 28	compressed-air driven homogenizing system	60418-00
	or	
XP44	compressed-air driven homogenizing system Instruction Manual	60417-00

### 2.2 Homogenising tools and accessories (to be ordered separately):

Tools	Description	Part No.
T 6	Shaft, 6 mm diameter, 55 mm long, complete with generator	60420*
T 6	Shaft, 6 mm diameter, 120 mm long, complete with generator	60410*
T 10	Shaft, 10 mm diameter complete with generator	60421*
T 17	Shaft, 10 mm diameter, complete with 17 mm generator	60426*
T 20	Shaft, 20 mm diameter, complete with generator	60422*
G 20	Shaft, 20 mm diameter, complete with generator	60423*
T 30	Shaft, 20 mm diameter, complete with 30 mm generator	60424*
G 30	Shaft, 20 mm diameter, complete with 30 mm generator	60425*
Pedestal	platform dimensions:	60491-00

<b>Tools</b>	<b>Description</b>	<b>Part No.</b>
<b>Stand incl. rod</b>	330x200 mm, rod diameter: 16mm, length: 690mm	
<b>Cross over clamp</b>	clamp to attach the unit to a stand Material:	60492-00
<b>Cross over clamp</b>	clamp to attach the unit to a stand Material:	60492-00
<b>Universal key</b>	Universal rotor wrench for shafts 6 - 30 mm	60470-00
<b>Socket wrench</b>	for all shafts	60471-00
<b>Rotor wrench</b>	for shaft T 10	60472-00
<b>Rotor wrench</b>	for shaft T 20/G 20	60473-00
<b>Rotor wrench</b>	for shaft T 30/G 30	60476-00
<b>Rotor wrench</b>	for shaft T 40/G 40	60477-00

### **2.3 Recommended Accerssories not available at CAT**

Coupling with minimum ID 4.5 mm/3/16"  
1/2" Service unit (Capacity 5 - 8 micron)

**Important Note:**

After reading and understanding the instruction manual you may now start operating the unit.

**Important Note:**

Read this instruction manual carefully before operating the instrument. Should there be any additional questions, after reading these instructions, concerning the set-up, operation or warranty, please contact either your distributor, or the manufacturer at the following address

# CAT

Ingenieurbüro CAT  
M. Zipperer GmbH  
Etzenbach 16  
D-79219 Staufen  
+49-7636-7803-0

**Important Note:**

Please store the instruction manual in a place easily accessible to every user.

**Important Note:**

Make sure that the unit is to be operated under the conditions corresponding to that as described under the section titled "Technical Data"

## 3 Safety Information

### 3.1 Explanation of the safety information in the instruction manual

This symbol emphasizes the following notices, which have to be strictly observed by the user! Disregarding of warnings may result in impairment of serviceability as well as impairment of the user.

This symbol emphasizes prohibitions, which have to be strictly observed. Any disregard of the prohibitions may result in severe impairment of serviceability as well as severe impairment of the user.

This symbol emphasizes notices, which have to be observed by the user to assure a safe operation of the unit.

### 3.2 Explanation of the safety information on the device

**CAUTION:**

This symbol indicates to read the instruction manual carefully prior to operation of the instrument. Please mark points which require special attention in your field of application so they are not overlooked. Disregarding of warnings may result in impairment of serviceability as well as impairment of the user.

### 3.3 General Safety Information

**Attention:**

Please comply with all safety and accident-prevention regulations applicable to laboratory work.

**Attention:**

Only instructed users may operate the instrument.

**Important Note:**

Make sure that all safety information of the instrument respectively on the device itself are clearly visible during operation

**Attention:**

If noise level exceeds 85 dB (A) at the work station operator must use ear plug and/or ear mufflers.

**Attention:**

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

Ensure that the unit is firmly attached to a stand.


## 4 Set-up and Operation

### 4.1 Intended Use

The XP28/XP44 is a compressed-air driven homogenizing system which in connection with a homogenizing tool produces emulsions and dispersions.

### 4.2 Description XP44/XP28

### 4.3 Set-up

The compressed-air driven homogenizing system XP28/XP44 may only be connected to the compressed air ductwork system via a maintenance unit. Prior to connecting the XP28/XP44, please check whether the pressure of your maintenance unit is the same as that one stated on the rating plate. Secondly check whether the lubrication is guaranteed by enough oil in the maintenance unit.

#### 4.4 Air supply requirements

The DIN ISO quality specifications require that air is free of solid particles larger than 5 micron, remaining water content to be max. 6g/m<sup>3</sup> (pressure dew point +3°C) and air may contain 5mg/m<sup>3</sup> of oil.

Therefore please observe the following:


- Operating with dirty and unlubricated air reduces the systems life substantially. Before installing extension hoses, blow them clean and minimize, if at all necessary, coupling and keep them clean.
- Conserve  
Put some drops of oil in the air inlet and briefly blow through.

#### 4.5 Operating Conditions

The compressed-air driven homogenizing system XP28/XP44 is designed for an operating pressure of 6 bar at full load. Pressure less than this will reduce output and speed, while appreciably higher pressure will lead to premature wear and excessive speed.

Clean, water-free and lubricated air is essential for good performance and a long machine life.

#### 4.6 Assembly of support rod

Screw the support rod (5) into the threaded bore (4) of the unit.

#### 4.7 Attaching the drive unit to a pedestal stand


To ensure safe operation the drive units are mounted to a pedestal stand using a cross over clamp. Vibrations may loosen the holding screws. Please check periodically whether all holding screws are securely tight.

- Attach cross over clamp to the stand rod.
- Now attach the drive unit to the stand by inserting the support rod (5) into the free opening of the cross over clamp.
- Tighten the knurled screw of the cross over clamp.

## 4.8 Inserting the homogenizing tool

The homogenizing tool comes always assembled. The connecting adapter (drive pin (6)) is located on the top of the shaft and is used to attach the shaft to the drive motor. When inserting the homogenising tool, make sure that the coupling fits well. This can be easily detected by the "click" that is heard when the coupling shaft has reached the end of the coupling slot of the motor. Do not forget to securely fasten the holding screw (2).

Never run homogenising tools dry as the bearings and gaskets will be damaged if the generator is not cooled by the medium. The distance between the homogenising tool and the vessel bottom should not be less than 10mm. The filling of the medium may not be less than about 55 mm. To avoid a vortex it is necessary to insert the dispersing shaft out of the center of the vessel.

If these conditions are met the unit is ready for operation.

## 4.9 Operation XP44

Prior to connecting the drive unit to your maintenance unit check whether the pressure of the maintenance unit corresponds to the pressure indicated on the rating plate of the compressed-air driven XP44. Check whether the lubrication is guaranteed by enough oil in the maintenance unit.

Before starting the unit set the speed regulator to the lowest speed.

The drive unit is turned on and off by means of a rotary switch on the motor housing. Turn the upper part to the right.

Switch position > 0 - 44,000: Motor running

Switch position 0: Motor off.

The speed can be steplessly adjusted at the sliding switch

The maximum speed heavily depends on the media and the size of homogenising tools used since the speed cannot be automatically adjusted with this drive unit as with an electric motor.

**Please find following some typical values using XP44 with a 20mm tool.**

P/(bar)	1	2	3	4	5	6	7
n/1min	8000	15000	20000	25000	29000	32000	35000

#### 4.10 Operation XP28

Prior to connecting the drive unit to your maintenance unit check whether the pressure of the maintenance unit corresponds to the pressure indicated on the rating plate of the compressed-air driven XP28. Check whether the lubrication is guaranteed by enough oil in the maintenance unit.

Before starting the unit set the speed regulator to the lowest speed.

The drive unit is turned on and off by means of a rotary switch on the motor housing. Turn the upper part to the right.

Switch position > 0 - 28,000: Motor running

Switch position 0: Motor off.

The speed can be steplessly adjusted at the sliding switch

The maximum speed heavily depends on the media and the size of homogenising tools used since the speed cannot be automatically adjusted with this drive unit as with an electric motor.

**Please find following some typical values using XP28 with a 20mm tool.**

P/(bar)	1	2	3	4	5	6	7
n/1min	6000	10000	13000	15500	17000	19000	20000

#### 4.11 Exchanging the homogenizing shafts

Prior to exchange the homogenizing shaft disconnect the unit from the compressed air supply.

Loosen the holding screw and firmly pull out the shaft by hand. It will disconnect easily. When inserting another shaft, make sure that the coupling fits well. This can be easily detected by the "click" that is heard when the coupling shaft has reached the end of the coupling slot of the motor. Do not forget to securely fasten the holding screw (2).

#### 4.12 Exchanging the generator:

To exchange the generator (rotor and stator), always take the shaft out of the motor drive unit as described above. For undoing the rotor, use the two special keys (to be ordered separately). The socket wrench is used to counter-hold the shaft, and the rotor wrench to undo the rotor. To exchange the stator, remember that the fastening thread is counter-clockwise. You will find further detailed descriptions in the instruction manual which comes with each homogenising tool.

#### **Important Note**

Always unscrew the rotor first before unscrewing the stator.

## 5 Maintenance and Cleaning

### 5.1 Cleaning the drive unit XP44/XP28

To clean the drive use only water with a detergent that contains tensides or use isopropylalcohol for stubborn soiling.

### 5.2 Cleaning the homogenising tool

To avoid clogging clean shaft, generator and in case of G-shaft the sealing parts after each use. This is done by operating it in a solvent which dissolves substance residues and is not harmful to the gasket. This is usually sufficient to clean the generator.

**Chemical sterilization** may be also a method. General-purpose disinfectants such as formalin, alcohol, etc. may be used. It is important to remove disinfectant residues with sterilized water.

**Important Note:**

Make sure that the bearings, O-rings and gaskets are resistant to solvents.

**Sterilization with moist heat:**

This denotes the use of a steam jet pressurized to 2 bar at 120°C.

### 5.3 Maintenance of drive unit XP44/XP28

Check idle speed periodically

Even if the tool still gives satisfactory service have it checked by trained service technicians after 300 to 400 operating hours, or at least once a year.

**Attention:**

Do not open the instrument. Repairs are only to be carried out by trained service technicians.

### 5.4 Maintenance homogenizing tools

The gaskets in the homogenizing tools must be constantly monitored. In the event of leakage the suction effect of the rotating shaft can cause the medium to penetrate as far as the drive unit. If liquid emerges from the side hole at the top of the shaft tube stop work immediately and check the gaskets. The function of the generators depends on the condition of the sharp edges on the rotor and stator. These edges may be blunted very quickly in abrasive media, reducing the effectiveness of homogenizing.

## 6 Dismantling, Transport and Storage

### 6.1 Dismantling

1. Switch the unit off.
2. Disconnect the unit from the maintenance unit and the compressed air ductwork system.
3. Remove any glass beakers and any other equipment around the drive unit.
4. Remove the homogenizing tool by loosening the knurled screw (2).
5. Loosen the cross over clamp and remove the clamp from the support rod.
6. Unscrew the support rod.
7. Now you may remove the instrument from the working area.

### 6.2 Transport and Storage

#### Prior to transport:

Switch the instrument off and proceed with dismantling as described under “Dismantling”.

Place the instrument and its parts in its original packaging or another suitable container to protect it during transport. Close the packaging with adhesive tape.

Store the instrument in a dry environment. Please observe the specified conditions of the ambient (temperature and humidity).

Do not subject the instrument to mechanical shocks or vibration during transporting it.

In case you do not use the original packaging please mark the box with the following notes:

- Glass symbol (handle with care, fragile)
- Umbrella (keep dry)
- Content (list of content)
- Storage ambient:
  - Max. ambient temperature : RT to +40°C
  - Max. humidity: 80%

## 7 Disposal

Please dispose of used instruments and defective components at your local recycling collection point. Prior to disposal, sort according to materials: metal, glass, plastic, etc. Also be sure to dispose of the packing material in an environmental-friendly manner.

## 8 Warranty and Liability

The manufacturer agrees to either repair, or replace, at the manufacturer's discretion, any defects in materials or workmanship which develop within 24 months of the delivery of this product to the original user. In the event of replacement, the replacement unit will be guaranteed for the remainder of the original twenty-four (24) month period or ninety (90) days, whichever is longer.

If this product should require service, contact your local distributor or manufacturer for necessary instructions.

This guarantee will not apply if the defect or malfunction was caused by accident, neglect, unreasonable use or fitness for a particular purpose, which extend beyond the description and period set forth herein.

The manufacturer's sole obligation under this guarantee is limited to the repair or replacement of a defective product and the manufacturer shall not, in any event, be liable for any incidental or consequential damages of any kind, resulting from use or possession of the product.

### Attention:

The user has to determine, if the instrument is suitable for his specific application. If there are any further queries, contact your local dealer or the manufacturer direct.

## 9 Technical Data

### 9.1 Technical Data XP44

Type: XP44	Specifications
Operating Pressure*	6 bar (90 PSI)
No load speed (idle speed)	44,000 rpm
Power	250 Watt
Torque	0.22Nm
Air consumption at full load	0.38 m <sup>3</sup> /min
Air consumption at no load	0.52 m <sup>3</sup> /min
Inside diameter of hose	8 mm
Permissible ambient temperature	5 - 40 °C
Permissible humidity	80 % RH
Sound pressure level ISO 15744	72 dB(A)
Vibration ISO 8662	< 2.5 m/s <sup>2</sup>
Case dimensions	Length: 200 mm, Diameter: 35 mm
Weight without hose	0.710 kg

## 9.2 Technical Data XP28

Type: XP28	Specifications
Operating Pressure*	6 bar (90 PSI)
No load speed (Idle speed)	28,000 rpm
Power	370 Watt
Air consumption at full load	0.6 m <sup>3</sup> /min
Air consumption at no load	0.6 m <sup>3</sup> /min
Inside diameter of hose	8 mm
Permissible ambient temperature	5 - 40 °C
Permissible humidity	80 % RH
Sound pressure level ISO 15744	80 dB(A)
Vibration ISO 8662	< 2.5 m/s <sup>2</sup>
Case dimensions	Length: 230 mm, Diameter: 40 mm
Weight without hose	1.0 kg

## 10 Repairs

When returning instruments for repair that have come into contact with hazardous substances, please:

- Fill in attached "Repair Return Form"
- Provide precise information on the relevant medium
- Take protective measures to ensure the safety of our receiving and maintenance personnel
- Mark the package as appropriate for hazardous materials.

## 10.1 Repair Return Form

### CONTACT/USER INFORMATION

Contact:	Phone No.
Fax No.	Email:
Billing:	Shipping:
Company	Company
Address	Address

### INSTRUMENT INFORMATION

Model	Serial No.	
Please describe all problems/malfunctions		
Operating Conditions (please fill in if applicable)		
Ambient Temp.	Humidity	Speed
Load	Volume	Viscosity
Temperature in °C	Sample Temperature	Operating Time
Sample Description*		

\*NOTE: If the instrument was exposed to hazardous material, it must be decontaminated BEFORE returning it to CAT and an MSDS for hazardous material must be included with the instrument.

### RETURN SHIPPING

UPS <input type="checkbox"/>	Air Parcel Post <input type="checkbox"/>	Collect** <input type="checkbox"/>	Other** <input type="checkbox"/>
------------------------------	--	------------------------------------	----------------------------------


\*\*Your account number is required for UPS collect respectively the address and contact of your preferred forwarder if you choose any other transport means.

### PACKAGING INSTRUCTIONS TO RETURN A CAT INSTRUMENT FOR REPAIR

- ✓ Remove all accessories (e.g. homogniser tools, stirring paddles) from the instrument
- ✓ Clean excess testing material off the instrument/accessory
- ✓ Include MSDS sheets for all hazardous materials used with this instrument
- ✓ Pack the instrument in its original box. If the box is not available, take care to wrap the instrument and accessories with enough material to support them.
- ✓ DO NOT send pedestal stand unless there is a problem with the upright rod, clamp or base. If there is a problem with the stand remove the upright rod from the base and individually wrap each item to avoid contact with the instrument. (Applicable for overhead stirrers and homogenisers.)
- ✓ Pack the instrument and related items in a strong box for shipping. Mark the outside of the box with handling instructions.

Example: "Handle with care" or "Fragile- Delicate Instrument" and send to:

**Ingenieurbüro CAT, M. Zipperer GmbH, Service Department, Etzenbach 16, D-79219 Staufen**

Ingenieurbüro  
M. Zipperer GmbH

Etzenbach 16, D-79219 Staufen  
Tel.: +49 (0) 7636 78030  
Fax: +49 (0) 7636 780345  
E-Mail: [info@cat-ing.de](mailto:info@cat-ing.de)  
Internet: <http://www.cat-ing.de>

 Declaration of conformity

**We declare under our sole responsibility that this product (see label) corresponds to the regulations and conforms to the standards or standardized documents**

**2006/42/EC,**  
DIN EN ISO 12100-1:2004-04  
DIN EN ISO 12100-2:2004-04

**In the case of a modification of the unit which has not been agreed on with us, this declaration becomes null and void and the warranty expires.**

<b>Compressed air-driven unit XP44</b>
<b>Compressed air-driven unit XP28</b>

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